

NEUROLAB COMPANY PRODUCT LIST

<http://www.neurolab.ru/en>



Product code	KPF-01	KPF-01b	KPF-01c-2	KPF-01c
Product name (Russian)	БиоМышь Профессиональная	БиоМышь Исследовательская		БиоМышь Индивидуальная
Product name (English)	BioMouse Professional	BioMouse Research		BioMouse Personal Edition
Development year	2004	2006	2007	2007
Hardware				
Short description	Professional, with sensors built into mouse (or with external sensors) attached to the separate case	Double-PCB. Built into the mouse case. Transitional model.	Single-PCB. Combined with mouse.	Single-PCB. Combined with mouse.
Registered signals	PPG, GSR, ESR	PPG	PPG	PPG
Hardware driver	Driver for FTDI USB ↔ COM adaptor	Driver for FTDI USB ↔ COM adaptor	LibUSB / WinUSB	LibUSB / WinUSB
Software				
Version	Professional, 2.4.299			Home, 2.3.929
Supported languages	Russian only			Russian, English. More on request.*
Supported methods	19	15	15	1
System requirements				
CPU	x86-compatible or x86-64-compatible at 1000 MHz or higher.			
RAM	At least the value recommended for installed operating system.			
Disk drives	Up to 100 Mb of HDD space is required for installation, depending on the operating system's configuration. Additional space is required for results storage.			
Optical drive	Optical drive capable of reading CD-ROM disks.			
I/O ports	One free USB 1.1 or 2.0 port			
Display	Required: resolution 800x600, 16-bit color. Recommended: resolution 1024x768 or higher, 32-bit color.			
Supported operating systems	Windows 2000, XP (SP1, SP2), Windows 2003 Server, Windows Vista			Windows 2000 or higher (except server OSes)
Required application software	MS Access, MS Excel, version 2000 or higher			—
Price of one piece, warehouse in Moscow				Price of 1 pc in China
1 pc +	Contact us!	Contact us!	Contact us!	Contact us!
500 pcs +	Contact us!	Contact us!	Contact us!	Contact us!
1 000 pcs +	Contact us!	Contact us!	Contact us!	Contact us!
500 000 pcs +	Contact us!	Contact us!	Contact us!	Contact us!

* For details on obtaining BioMouse in languages other than Russian or English please contact us.